



Air Bedside Call Unit

ZKRUNCWRWL10100

ZKR

The Air Bedside Call Unit is used in patient rooms. There are call and cancel buttons on the unit. In an emergency, a patient uses the call button to make an emergency call which appears as an alert on the Nurse Panel. Typical locations for this unit are on the walls of patient rooms and living areas.

- Air Patient Handsets work through connection to this unit
- Can be flush mounted and has an ergonomic design
- Battery lasts up to 5000 calls
- Made of high quality PC-ABS-FR Raw Materials



Technical Specification

ZKRUNCWRWL10100

<i>Weight</i>	60 g
<i>Width</i>	80 mm
<i>Depth</i>	35 mm
<i>Height</i>	80 mm
<i>Mounting Dimensions</i>	W-55mm D-20mm H-55mm
<i>Processor</i>	ARM
<i>Power Supply</i>	12V Battery
<i>Operating Temperature / Humidity</i>	-10°C +60°C / % 30 %80 RH
<i>Raw Material</i>	PC-ABS-FR
<i>Handset Connection</i>	1xRJ45
<i>Battery Life</i>	5000 calls
<i>Low Battery Alert</i>	Yes
<i>Color</i>	White
<i>Wireless Range</i>	25 m
<i>Reassurance Light</i>	Yes
<i>Wireless Frequency</i>	433 Mhz Encrypted Communication Protocol