



Code Blue Unit

ZKRUNCIPCDBL101

ZKR

The Code Blue Unit is mounted on the wall as needed in the hospital. There are backlit call and code blue buttons on the unit. The call button may be pressed to initiate a nurse call which appears as an alert on the Nurse Control Panel.

The code blue button may be pressed to start a Code Blue call. This button can also be configured as an assist call button.



Technical Specification

ZKRUNCIPCDBL101

<i>Weight</i>	56 g
<i>Width</i>	80 mm
<i>Depth</i>	35 mm
<i>Height</i>	80 mm
<i>Mounting Dimensions</i>	W-55mm D-20mm H-55mm
<i>Processor</i>	ARM
<i>Operating Voltage</i>	12V
<i>Operating Temperature / Humidity</i>	-20°C +60°C / % 30 %80 RH
<i>Raw Material</i>	PC-ABS-FR
<i>Buzzer Speaker</i>	2400Hz 65 dB
<i>Ports and Slots</i>	2xRJ45
<i>Color</i>	White
<i>Back Light</i>	LED