


Senset 3 Input Module

ZKRUNCIPIZLMOD1

ZKR

Provides interface with various external equipment such as ECG Patient Monitors or Saline Solution bags to track the quantity of remaining Saline Solution. This module allows the Nurse Call System to generate notifications and alarms automatically through data received from the existing medical equipment and devices. This way healthcare personnel are notified immediately even when the patient in need of assistance is unable to make a call manually.



Technical Specification

ZKRUNCIPIZLMOD1

Weight	100 g
Width	34 mm
Depth	48 mm
Height	20 mm
Operating Temperature	-20°C +60°C
Ports and Slots	2xRJ45, 3xDry Contact
Operating Voltage	12 V
Color	Black